

Appendix I

*Gynecological Papanicolaou Staining Procedure

Stain Station	Reagent	Suggested Time
1	95% Ethanol	3 min
2	70% Ethanol	1 min
3	50% Ethanol	1 min
4	Distilled Water	1 min
5	Hematoxylin I	7 min
6	Distilled Water	10 sec
7	Clarifier I	1 min
8	Distilled Water	1 min
9	Bluing Solution	30 sec
10	Distilled Water	30 sec
11	50% Ethanol	30 sec
12	95% Ethanol	30 sec
13	OG-6 Solution	2 min
14	95% Ethanol	15 sec
15	95% Ethanol	15 sec
16	EA-50 Solution	4 min
17	95% Ethanol	1 min
18	95% Ethanol	1 min
19	100% Ethanol	30 sec
20	100% Ethanol	30 sec
21	100% Ethanol	30 sec
22	Xylene	1 min
23	Xylene	3 min
24	Xylene	1 min

*Please note: The Staining Procedures in Appendix I & II are only suggested staining procedures.

chemQ bioscience LLC
Research Triangle Park, NC 27709 USA
techsupport@vitroprep.com

Appendix II

*Non-Gynecological Papanicolaou Staining Procedure

Stain Station	Reagent	Suggested Time
1	95% Ethanol	3 min
2	95% Ethanol	3 min
3	Distilled Water	30 sec
4	Hematoxylin 7211	40 sec
5	Distilled Water	90 sec
6	Clarifier I	30 sec
7	Distilled Water	1 min
8	Bluing Solution	1 min
9	Distilled Water	1 min
10	95% Ethanol	30 sec
11	OG-6 Solution	2 min
12	95% Ethanol	15 sec
13	95% Ethanol	15 sec
14	EA-50 Solution	3 min
15	95% Ethanol	15 sec
16	95% Ethanol	15 sec
17	100% Ethanol	1 min
18	100% Ethanol	1 min
19	100% Ethanol	1 min
20	Xylene	1 min
21	Xylene	1 min
22	Xylene	1 min

*Please note: The Staining Procedures in Appendix I & II are only suggested staining procedures.

chemQ bioscience LLC
Research Triangle Park, NC 27709 USA
techsupport@vitroprep.com