

## ***VitroPrep™ Cytology Processing Kit***

- ✓ Prepares clinical specimens for unobstructed clear analysis of cells by pathologists.
- ✓ Preserves the cells, lysis red blood cells & removes mucus and debris.
- ✓ Provides an unobstructed mono-layer cell adhesion to standard microscope slides.
- ✓ Formulated using “**Green Chemistry**” technology.
- ✓ Utilizes a simple processing procedure that requires no automation.
- ✓ Product of USA, CE certified, and ISO 13485 compliant.
- ✓ **Cost effective and affordable.**

❖ ***The Kit consists of two proprietary solutions***

### **PreservPlus™ & CytoBase™**



29 PreservPlus™ & 1 CytoBase™ Solutions

## ***PreservPlus™ Solution***



- ✓ A proprietary solution formulated to maintain cell morphology and DNA/RNA/protein quality.
- ✓ A new, gentle cleaning and RBC lysing solution in a single bottle.
- ✓ Formulated using “**Green Chemistry**” technology.

## ***CytoBase™ Solution***



- ✓ A proprietary solution which provides mono-layer cell adhesion to standard glass slides.
- ✓ No heating is required after treatment of the slides by the CytoBase™ Solution.
- ✓ Compatible with many standard stains.
- ✓ Formulated using “**Green Chemistry**” technology.

## ***VitroPrep™ Cytology Processing Procedure***

- ✓ Utilizes a simple processing procedure that can be adapted to any laboratory worldwide.
- ✓ Requires no automated instrumentations.
- ✓ Uses only a vortex and a standard swing-bucket centrifuge.


### *Follow instructions below for sample processing:*

1. Vortex PreservPlus™ Solution containing specimen\* for about 10 seconds.
2. Pour the specimen into a 15 mL conical centrifuge tube.
3. Centrifuge the tube for 10 min @ 1000 ×g.
4. Gently and quickly decant supernatant liquid from the tube.
5. Depending on the cellular pellet size, add CytoBase™ Solution in a ratio of 3 parts CytoBase™ Solution to 1 part of cellular pellet, typically 250-300 µL, into the tube and vortex the tube about 20 seconds.
6. Withdraw 45 µL of the mixture (step 5) and spread gently on a standard glass slide.
7. Allow drying of the slide at room temperature for 3-4 hours then follow standard Papanicolaou staining protocol.


\*For specimens containing blood, allow PreservPlus™ Solution containing specimen to sit for 24 hours for effective hemolysis.

## ***VitroPrep™ Clinical Slide Images***


*Below are some examples of study slides using VitroPrep™ Cytology Processing Kit.*



Negative for intraepithelial lesion & malignancy



Low grade squamous intraepithelial lesion



High grade squamous intraepithelial lesion